

The International Brain Laboratory

Brain-Wide Map Data – Hands-On Training

295501 neurons
547 insertions
194 brain regions
115 mice performing the task

Karolina Socha
PostDoc in Carandini and Harris Labs




INTERNATIONAL
BRAIN
LABORATORY

NeuroDataShare 2023
22 February 2023, UCL-SWC, London

The IBL: 22 labs working to understand decision-making

INTERNATIONAL
BRAIN
LABORATORY




 **FLATIRON**
INSTITUTE
Division of Simons Foundation

SIMONS
FOUNDATION

wellcometrust

1. **Brain-Wide Map (BWM) Data** (recordings from 194 brain regions during perceptual task)
2. Reproducible Ephys Data (recordings in the same location)
3. Behavior Data (only behavior during learning the task)
4. Future: individual projects (mesoscale imaging, photometry, optogenetics)

Future Tutorials: COSYNE 2023



INTERNATIONAL
BRAIN
LABORATORY

DATE: 9 March 2023

TIME: 10:00 - 11:00

LOCATION: Av. Duluth, Convention Floor, Fairmont The Queen Elizabeth

Access to data:

<https://viz.internationalbrainlab.org/>

Methods:

[https://figshare.com/articles/preprint/Data release
- Brainwide map - Q4 2022/21400815](https://figshare.com/articles/preprint/Data_release_-_Brainwide_map_-_Q4_2022/21400815)

Data Type:

[https://docs.google.com/document/d/1OqlqqakPak
HXRAwceYLwFY9gOrm8_P62XIfCTnHwstg/edit#](https://docs.google.com/document/d/1OqlqqakPakHXRAwceYLwFY9gOrm8_P62XIfCTnHwstg/edit#)


Example scripts:

[https://github.com/int-brain-
lab/UCL_NeuroDataShare2023](https://github.com/int-brainlab/UCL_NeuroDataShare2023)

1. **Experimental design**
2. **Data collection strategy**
3. **Data Integration**
4. **Getting started:**
installation
loading data
searching sessions
5. **Summary**

Experimental setup


3 cameras allowing to record behavioral responses
(licking, face motion, movement)




eid: 1a166c4f-4c53-422e-8473-b3cff85e6750
danlab/DY_009/2020-02-29/001

IBL decision-making task


Visual perceptual task



50-50 block (first 90 trials)




20-80 block (from 90th trial)



Task allows to measure neural correlates of:

- Sensory stimuli (visual detection)
- Motor response (wheel movement, or direction of movement)
- Reward signal
- Bias/Prior (at low stimuli contrasts)

IBL passive (after the task, data not released)




Protocol allows to measure:

- Post-task activity (spontaneous)
- Receptive fields mapping (visual responses)
- Sensory task-related responses (valve clicks, sound cues, visual stimuli)
- Spontaneous movement, licking
- Brain state modulations (activity during task vs passive)


Recording neural activity during behavior

Neuropixels probes
385 active channels




Recordings strategy



Targeted insertion
coordinates



Targeted recording sites



Example recordings
sites



Recordings acquisition:
spikeGLX (Bill Karsh, HHMI)

BWM Data Release:


295501 neurons

547 insertions

194 brain regions

115 mice

Data flow within and outside IBL



1. Raw data

- Subject-level metadata (Alyx database)
- Behavioral performance (*.json)
- Electrophysiology recordings (compressed *.cbin, *.meta)
- Behavioral video (compressed *.mp4)
- Histological probe track reconstruction (Alyx database)
- Ambient data (temperature, humidity, etc)

2. Post-processing data (extraction spikes, task signals, and video details)

- Spike sorting (pykilosort)
- Video tracking (DeepLabCut)
- Task signals extraction (wheel, photodiode, task event signals)
- Synchronization with the probe

3. Quality Control

- Hardware
- Electrophysiology recording, videos
- Behavioral performance
- Single units

IBL folder structure

lab / mouse / day / session / folder_name


cortexlab / KS091 / 2022-07-06 / 001 / folder_name

folder_name can be:

- ***alf*** – extracted data, including spike-sorted data
- *raw_ephys_data/probe00* ,
raw_ephys_data/probe01
- *raw_video_data*

```
cortexlab/KS091
|--2022-07-06/
|   |-- 001/
|       |--alf/
|           |-- probe00/
|           |-- pykilosort/
|           |-- probe01/
|           |-- pykilosort/
|           |-- raw_ephys_data/
|           |-- probe00/
|           |-- probe01/
|           |-- raw_video_data
```

IBL data type: ALF (ALyx Files)



Filename	Size
<code>spikes.times.npy</code>	<code>nSpikes</code>
<code>spikes.amps.npy</code>	<code>nSpikes</code>
<code>spikes.clusters.npy</code>	<code>nSpikes</code>
<code>clusters.amps.npy</code>	<code>nClusters</code>
<code>clusters.mlappv.npy</code>	<code>nClusters x 3</code>

Standard attributes

Event series	<code>*.times *.* _times[_timescale]</code> e.g. <code>trials.reward_times</code>
Interval series	<code>*.intervals[_timescale]</code> e.g. <code>tones.intervals</code>
Continuous timeseries	<code>*.timestamps</code> e.g. <code>video_trialStart.timestamps</code>

IBL data objects: *trials*, *wheel*, *camera*, *spikes*, *clusters*


<https://docs.google.com/document/d/1OqIqqakPakHXRAwceYLwFY9gOrm8P62XIfCTnHwstg/edit#heading=h.nvzaz0fozs8h>

NPY in matlab: <https://github.com/kwikteam/npy-matlab/>


1. **Raw data (compressed; bin files (mtscomp)**
<https://github.com/int-brain-lab/mtscomp>; avi video (ffmpeg))
2. **ALF (ALyx Files, extracted data)**
 - A. Post-processing session data, specific to session**
 - Trial data: trials.table
 - Wheel data: wheel.position, wheel.timestamps
 - DLC data: camera.dlc, camera.times, camera.ROIMotionEnergy, ROIMotionEnergy.position, licks.times
 - Passive data: _ibl_passiveRFM.times, _ibl_passiveStims.table, _ibl_passiveGabor.table
 - B. Post-processing ephys data, specific to probe**
 - Ephys data: _spikeglx_sync.channels, _spikeglx_sync.times, _iblqc_ephysSpectralDensity.power, _iblqc_ephysTimeRms.rms
 - Spike-sorted data: spikes.times, spikes.clusters, spikes.depths, clusters.metrics, clusters.acronym

Task trial structure and wheel data

Trial structure



Example neural activity during a single trial



trials.table:


feedback_times,
firstMovement_times, *stimOn_times*,
goCue_times

wheel:


_ibl_wheel.position : absolute linear displacement of wheel (cm).
_ibl_wheel.timestamps

Task data – single trial details

Example neural activity during a single trial



Single unit activity aligned to task events



Cluster #493, Caudoputamen

timestamps:

feedback_times, *firstMovement_times*, *stimOn_times*, *goCue_times*

details :

choice : -1 (turn **CCW**), +1 (turn **CW**), or 0 (**no-go**)

contrastLeft : **contrast** of left-side stimulus **0**, **6.25%**, **12.5%**, **25%**, **100%**


contrastRight : contrast of right-side stimulus (nan if stimulus is on the other side)

feedbackType : -1 for negative (**incorrect**), 1 for positive (**correct**), 0 for **no-go**


probabilityLeft : **0.5**, **0.2**, **0.8**

DLC data from videos

Right camera (150Hz)
Left camera (60Hz)



Body camera (30 Hz)



DLC data


[Details](#)

Dataset type	Collection	Dataset name
camera.dlc - LEFT	alf	_ibl_leftCamera.dlc.pqt
camera.dlc - RIGHT	alf	_ibl_rightCamera.dlc.pqt
camera.dlc - BODY	alf	_ibl_bodyCamera.dlc.pqt
camera.times - LEFT	alf	_ibl_leftCamera.times.npy
camera.times - RIGHT	alf	_ibl_rightCamera.times.npy
camera.times - BODY	alf	_ibl_bodyCamera.times.npy
camera.ROIMotionEnergy - LEFT	alf	leftCamera.ROIMotionEnergy.npy
camera.ROIMotionEnergy - RIGHT	alf	rightCamera.ROIMotionEnergy.npy
camera.ROIMotionEnergy - BODY	alf	bodyCamera.ROIMotionEnergy.npy
ROIMotionEnergy.position - LEFT	alf	leftROIMotionEnergy.position.npy
ROIMotionEnergy.position - RIGHT	alf	rightROIMotionEnergy.position.npy
ROIMotionEnergy.position - BODY	alf	bodyROIMotionEnergy.position.npy
camera.features - LEFT	alf	_ibl_leftCamera.features.pqt
camera.features - RIGHT	alf	_ibl_rightCamera.features.pqt
licks.times	alf	licks.times.npy

Raw video data

_iblrig_Camera.raw - LEFT	raw_video_data	_iblrig_leftCamera.raw.mp4
_iblrig_Camera.frameData - LEFT	raw_video_data	_iblrig_leftCamera.frameData.bin

Example plots using DLC variables



How to find session, or specific probe ?

Every session and probe have a unique number:

Probe ID (pid): *00b05238-aa75-4846-a480-c5ffef4529dc*

- Brain region location: clusters.acronym
- Ephys data: spikes.times, spikes.clusters

Experiment ID (eid): *258b4a8b-28e3-4c18-9f86-1ea2bc0dc806*

- task trials data: trials.table
- wheel data: wheel.position, wheel.timestamps
- dlc data: camera.dlc, camera.times, licks.times, camera.ROIMotionEnergy
- passive data: _ibl_passiveRFM.times, _ibl_passiveStims.table, _ibl_passiveGabor.table

Getting started

- access public data with ONE

IBL promotes open science, hence developed tools and packages are in **python**



+




colab

Amazon Web
Service



OpenAlyx



bulk data

metadata data

Data Access

**ONE
protocol**



You

```
one.load_collection(eid, 'alf', download_only=True)
one.load_object(eid, 'trials', collection='alf')
one.load_dataset(eid, '_ibl_wheel.timestamps.npy')
```


Getting started

- example notebook scripts

1. Installation
2. Searching sessions
3. Downloading sessions
4. Loading spikes, and trials to jupyter notebook
5. Query session
6. Plotting data

Documentation:

https://int-brain-lab.github.io/iblenv/notebooks_external/data_download.html

Getting started

- creating env and installation

Step01: Create environment

```
conda create --name ibl python=3.9  
conda activate ibl
```

Step02: Install IBL packages

```
pip install ONE-api  
pip install ibllib
```

ONE-api: *a set of functions that allow you to search the database and download the bulk data to your local computer*

Step 03: Setting up credentials

```
from one.api import ONE  
pw = 'international'  
one = ONE(base_url='https://openalyx.internationalbrainlab.org',  
password=pw, silent=True)
```

Getting started

- cloning BWM github repository

Step04: Clone BWM git repository

```
git clone https://github.com/int-brain-lab/paper-brain-wide-map.git
cd paper-brain-wide-map
pip install -e .
```

List of all sessions:

*paper-brain-wide-map / brainwidemap / fixtures / **2022_10_bwm_release.csv***

List of terms: ['pid', 'eid', 'probe_name', 'session_number', 'date', 'subject', 'lab']

	pid	eid	probe_name	session_number	date	subject	lab
0	56f2a378-78d2-4132-b3c8-8c1ba82be598	6713a4a7-faed-4df2-acab-ee4e63326f8d	probe00	1	2020-02-18	NYU-11	angelakilab
1	47be9ae4-290f-46ab-b047-952bc3a1a509	56956777-dca5-468c-87cb-78150432cc57	probe01	1	2020-02-21	NYU-11	angelakilab
2	6be21156-33b0-4f70-9a0f-65b3e3cd6d4a	56956777-dca5-468c-87cb-78150432cc57	probe00	1	2020-02-21	NYU-11	angelakilab
...
542	8bf0f1a4-0d8c-4df3-a99e-f7c81c809652	993c7024-0abc-4028-ad30-d397ad55b084	probe01	1	2020-09-16	CSH_ZAD_029	zadorlab

Getting started

- cloning BWM github repository

Step05: Clone git repository with examples for
UCL_NeuroDataShare2023

```
git clone https://github.com/int-brain-lab/UCL_NeuroDataShare2023.git
```

Example01_LaunchONE-api_SearchSe...

Example02_DownloadSession.ipynb

Example03_DownloadRawData.ipynb

Example04_LoadTrials_Spikes.ipynb

Example05_QuerySessions_BrainRegi...

Example06_QueryNeurons_GoodIBLU...

Example07_LoadWheel_plotTrace.ipy...

Example08_LoadData_PlotPSTH.ipynb

Example09_LoadData_PlotFR.ipynb

Example10_ReactionTimes.ipynb

Mayo Faulkner (IBL-BWM tutorials)

Example01: Launch ONE-api and search ephys sessions

1. Using file **2022_10_bwm_release**

```
from brainwidemap import bwm_query
import numpy as np
from one.api import ONE
one = ONE()
ba = AllenAtlas()

# load info about all released sessions
bwm_df=bwm_query(one=None, alignment_resolved=True, return_details=False, freeze='2022_10_bwm_release')
```

Loading bwm_query results from fixtures/2022_10_bwm_release.csv

bwm_df									
		pid		eid	probe_name	session_number	date	subject	lab
0	56f2a378-78d2-4132-b3c8-8c1ba82be598	6713a4a7-faed-4df2-acab-ee4e63326f8d			probe00	1	2020-02-18	NYU-11	angelakilab
1	47be9ae4-290f-46ab-b047-952bc3a1a509	56956777-dca5-468c-87cb-78150432cc57			probe01	1	2020-02-21	NYU-11	angelakilab
2	6be21156-33b0-4f70-9a0f-65b3e3cd6d4a	56956777-dca5-468c-87cb-78150432cc57			probe00	1	2020-02-21	NYU-11	angelakilab
3	1e176f17-d00f-49bb-87ff-26d237b525f1	a8a8af78-16de-4841-ab07-fde4b5281a03			probe00	1	2020-01-22	NYU-12	angelakilab
4	701026df-e170-4ca7-88aa-eb0b95ef6ba1	a8a8af78-16de-4841-ab07-fde4b5281a03			probe01	1	2020-01-22	NYU-12	angelakilab
...
542	8bf0f1a4-0d8c-4df3-a99e-f7c81c809652	993c7024-0abc-4028-ad30-d397ad55b084			probe01	1	2020-09-16	CSH_ZAD_029	zadorlab
543	5d570bf6-a4c6-4bf1-a14b-2c878c84ef0e	fece187f-b47f-4870-a1d6-619afe942a7d			probe01	1	2020-09-17	CSH_ZAD_029	zadorlab
544	f7c93877-ec05-4091-a003-e69fae0f2fa8	fece187f-b47f-4870-a1d6-619afe942a7d			probe00	1	2020-09-17	CSH_ZAD_029	zadorlab
545	675952a4-e8b3-4e82-a179-cc970d5a8b01	c7bd79c9-c47e-4ea5-aea3-74dda991b48e			probe01	1	2020-09-19	CSH_ZAD_029	zadorlab
546	79f44ba1-c931-4346-82eb-f628a9374045	c7bd79c9-c47e-4ea5-aea3-74dda991b48e			probe00	1	2020-09-19	CSH_ZAD_029	zadorlab

547 rows x 7 columns

Example01: Launch ONE-api and search ephys sessions

2. Using **one.search**

```
example_sess.keys()
```

```
dict_keys(['id', 'subject', 'start_time', 'number', 'lab', 'projects', 'url', 'task_protocol'])
```

subject = 'SWC_054'

sessions = **one.search**(subject=subject)

- this gives a list of only **eids** of all detected sessions

```
In [58]: ## searching sessions from sepcific subject  
         subject = 'SWC_054'  
         # query sessions endpoint  
         sessions = one.search(subject=subject)  
         print(f'No. of detected sessions: {len(sessions)}'))  
         sessions
```

```
No. of detected sessions: 57
```

```
Out[58]: ['25731502-95bd-4aa7-b5e9-87414a3c4be6',  
          '6bb5da8f-6858-4fdd-96d9-c34b3b841593',  
          '671c7ea7-6726-4fbc-adeb-f89c2c8e489b',  
          'eebacd5a-7dcd-4ba6-9dff-ec2a4d2f19e0',  
          '5c7d2345-1f0e-40e5-aad7-2c6133b71b09',  
          '6c6983ef-7383-4989-9183-32b1a300d17a',
```


INTERNATIONAL
BRAIN
LABORATORY

```
insertions = one.alyx.rest('insertions', 'list', subject=subject)
```

- this gives **metadata** about session

```
Out[59]: [{'id': '7909c0aa-c074-4e19-aabf-b8167c682a5b',
'session': '6bb5da8f-6858-4fdd-96d9-c34b3b841593',
'model': '3B2',
'session_info': {'subject': 'SWC_054',
'start_time': '2020-10-11T20:00:42.054571',
'number': 1,
'lab': 'mrsicflogellab',
'id': '6bb5da8f-6858-4fdd-96d9-c34b3b841593',
'task_protocol': '_iblrig_tasks_ephysChoiceWorld6.4.2'},
'name': 'probe00',
'json': {'qc': 'WARNING',
'n_units': 139,
'xyz_picks': [[-2013, -3024, -43],
[-2013, -3099, -292],
[-2013, -3125, -393],
[-2013, -3149, -568],
```

Example02: Navigate and load data

Use **bwm_query** to generate data frame table with IBL sessions list and get eids and pids from this table:

```
In [11]: from brainwidemap import bwm_query
import numpy as np
from one.api import ONE
one = ONE()
ba = AllenAtlas()

# load info about all released sessions
bwm_df=bwm_query(one=None, alignment_resolved=True, return_details=False, freeze='2022_10_bwm_release')

Loading bwm_query results from fixtures/2022_10_bwm_release.csv
```

```
In [10]: bwm_df
```

```
Out[10]:
```

	pid		eid	probe_name	session_number	date	subject	lab
0	56f2a378-78d2-4132-b3c8-8c1ba82be598	6713a4a7-faed-4df2-acab-ee4e63326f8d		probe00	1	2020-02-18	NYU-11	angelakilab
1	47be9ae4-290f-46ab-b047-952bc3a1a509	56956777-dca5-468c-87cb-78150432cc57		probe01	1	2020-02-21	NYU-11	angelakilab

Example02: Download data

1. Using `one.load_collection`

```
eid='288bfbf3-3700-4abe-b6e4-130b5c541e61'  
sessions = one.load_collection(eid, 'alf', download_only=True)
```

The format of the returned datasets gives the path of the collection followed by the dataset. e.g in the case of `alf/trials.table.pqt`, **alf** is the collection and `trials.table.pqt` is the dataset.

The collection is important as it differentiates datasets with the same name e.g `spikes.times` in `alf/probe00/pykilosort` and `spikes.times` in `alf/probe01/pykilosort`.

```
collections = one.list_collections(eid)  
print(collections)
```

```
['alf/probe01/pykilosort', 'alf/probe00', 'alf/probe00/pykilosort', 'alf/probe01', 'raw_ephys_data/probe00',  
obe01', 'alf', 'raw_passive_data', 'raw_ephys_data', 'raw_behavior_data', 'spike_sorters/pykilosort/probe01',  
'spike_sorters/pykilosort/probe00', 'spike_sorters/ks2_matlab/probe00', 'spike_sorters/ks2_matlab/probe01']
```

Example02: Download data

2. Using `one.load_dataset`

```
eid='288bfbf3-3700-4abe-b6e4-130b5c541e61'
```

```
spike_times = one.load_dataset (eid, 'spikes.times.npy',  
                                collection='alf/probe00/pykilosort')
```

A single dataset can be downloaded and loaded into memory by passing in the eid and dataset as arguments into the `one.load_dataset` method,

```
# Download and load the left camera timestamps  
left_cam_times = one.load_dataset(eid, '_ibl_leftCamera.times.npy')  
  
# Download and load the spikes times for probe00  
spike_times = one.load_dataset(eid, 'spikes.times.npy', collection='alf/probe00/pykilosort')
```

Example02: Download data

3. Using **one.load_object**

eid='288bfbf3-3700-4abe-b6e4-130b5c541e61'

trials = one.load_object(eid, 'trials', collection='alf')

A group of attributes (e.g amps, depths, metrics) belonging to the same object (e.g clusters) can be downloaded and loaded in one command using the `one.load_object` method

```
# Load in all trials datasets
trials = one.load_object(eid, 'trials', collection='alf')
wheel = one.load_object(eid, 'wheel', collection='alf')

# Only download the clusters object for probe01
clusters = one.load_object(eid, 'clusters', collection=f'alf/{pname}/pykilosort', download_only=True)

# Only download the spikes object for probe01
spikes = one.load_object(eid, 'spikes', collection=f'alf/{pname}/pykilosort', download_only=True)
```

Example03: Download raw data

Use **one.load_datasets** to load *lf.cbin , *ap.cbin data (use PID)

```
In [*]: from one.api import ONE
import spikeglx
one = ONE()

pid = 'da8dfec1-d265-44e8-84ce-6ae9c109b8bd'
eid, probe = one.pid2eid(pid)

band = 'ap' # either 'ap', 'lf'

# Find the relevant datasets and download them
dsets = one.list_datasets(eid, collection=f'raw_ephys_data/{probe}', filename='*.lf.*')
data_files, _ = one.load_datasets(eid, dsets, download_only=True)
bin_file = next(df for df in data_files if df.suffix == '.cbin')

K:\Flatiron\ONE\alyx.internationalbrainlab.org\hoferlab\Subjects\SWC_043\2020-09-21\001\raw_ephys_data\probe00\spikeglx
```

Use **one.load_datasets** to load *mp4 video files (use EID)

```
In [ ]: from one.api import ONE
import ibllib.io.video as vidio

one = ONE()
eid = '4ecb5d24-f5cc-402c-be28-9d0f7cb14b3a'
label = 'body' # 'left', 'right' or 'body'

# Load raw video
video_body = one.load_dataset(eid, f'{label}Camera.raw*', collection='raw_video_data')
```

More examples: https://int-brain-lab.github.io/iblenv/loading_examples.html

Example04: Load session data to jupyter notebook

Use **SessionLoader** to load e.g. trials, wheel data

```
# import session loader
from one.api import ONE
from brainbox.io.one import SessionLoader
one = ONE()
eid='6713a4a7-faed-4df2-acab-ee4e63326f8d'
# instantiate session loader
sess_loader = SessionLoader(one=one, eid=eid)
```

```
SessionLoader(one=One (online, https://alyx.internationalbrainlab.org), session_path=WindowsPath('K:/Flatiron/ONE/alyx.internationalbrainlab.org/angelakilab/Subjects/NYU-11/2020-02-18/001'), eid='6713a4a7-faed-4df2-acab-ee4e63326f8d')
```

```
# Load in trials data
sess_loader.load_trials()
```

```
# Load in wheel data
sess_loader.load_wheel()
```

```
sess_loader.data_info
```

```
sess_loader.trials.keys()
```

```
Index(['stimOff_times', 'goCueTrigger_times', 'intervals_bpod_0',
       'intervals_bpod_1', 'firstMovement_times', 'goCue_times',
       'probabilityLeft', 'response_times', 'feedbackType', 'rewardVolume',
       'contrastRight', 'choice', 'feedback_times', 'stimOn_times',
       'contrastLeft', 'intervals_0', 'intervals_1'],
      dtype='object')
```

Example04: Load session data to jupyter notebook

trials table structure:

sess_loader.trials

	stimOff_times	goCueTrigger_times	intervals_bpod_0	intervals_bpod_1	firstMovement_times	goCue_times	probabilityLeft	response_times	feedbackType
0	315.645189	268.879047	0.000000	62.662702	312.82632	268.879856	0.5	313.593986	-1.0
1	377.427764	317.356906	63.199299	126.375502	NaN	317.357932	0.5	377.358385	-1.0
2	386.672888	381.061309	126.884199	133.690002	385.42132	381.062242	0.5	385.593944	1.0
3	406.660838	389.144464	134.121899	153.677802	404.07032	389.145186	0.5	404.595102	-1.0

Example04: Load session data to jupyter notebook

Use **SpikeSortingLoader** to load spikes data

The SpikeSorting loader can be used in spike sorting data for a single insertion. It can be instantiated with an ONE instance and either a **pid** or **and eid, pname** combination

```
pid='56f2a378-78d2-4132-b3c8-8c1ba82be598'
pname='probe00'
# instantiate with a pid
spike_loader = SpikeSortingLoader(pid=pid, one=one)

# alternatively instantiate with an eid and probe name
spike_loader = SpikeSortingLoader(eid=eid, one=one, pname=pname)

# Download and load data
spikes, clusters, channels = spike_loader.load_spike_sorting()
# Assign brain location information from channels to clusters
clusters = spike_loader.merge_clusters(spikes, clusters, channels)
spikes

{'depths': array([439.90316772, 273.93023682, 445.34341431, ..., 28.31001663,
191.40383911, 41.63357162]),
'clusters': array([ 99,  64, 324, ...,   6,  28,   9], dtype=uint32),
'times': array([8.38230849e-03, 1.27823046e-02, 1.41823034e-02, ...,
4.68878040e+03, 4.68878233e+03, 4.68878343e+03]),
'amps': array([1.48818473e-04, 7.21604797e-05, 7.99972536e-05, ...,
1.45808202e-04, 8.26257018e-05, 1.57378919e-04])}

clusters.keys()

dict_keys(['uids', 'depths', 'channels', 'cluster_id', 'amp_max', 'amp_min', 'amp_median', 'amp_std_db', 'contamination', 'contamination_alt', 'drift', 'missed_spikes_est', 'noise_cutoff', 'presence_ratio', 'presence_ratio_std', 'slidingRP_viol', 'spike_count', 'firing_rate', 'label', 'x', 'y', 'z', 'acronym', 'atlas_id', 'axial_um', 'lateral_um'])
```

Example05: Query sessions

Use **one.alyx.rest** to query sessions from specific brain region

```
In [16]: # Find sessions recorded in specific brain region

# Loading data with SC
#from oneibl.one import ONE
from one.api import ONE
import matplotlib.pyplot as plt
import pandas as pd
import numpy as np
one = ONE()

## FIND EIDS IN SPECIFIC BRAIN REGIONS
# example brain regions
#brainregions_acronyms=["LGD", "CP", "MOp", "VISp", "ZI", "SNr"]
brain_region="LGD"
insertions = one.alyx.rest('insertions', 'list', task_protocol='ephys',
                           atlas_acronym=brain_region,
                           project='ibl_neuropixel_brainwide_01', no_cache=True)

probe_insertions = [p['name'] for p in insertions]
eid_insertions = [s['session'] for s in insertions]
pid_insertions=[p['id'] for p in insertions]
subject_insertions=[m['session_info']['subject'] for m in insertions]
start_time_insertions = [k['session_info']['start_time'][0:10] for k in insertions]

data_in={'subject':subject_insertions,'day':start_time_insertions,'probe':probe_insertions,
        'eid':eid_insertions,'pid':pid_insertions,'brain_region':brain_region}

# create data frame
df_experiments=pd.DataFrame(data_in)

print('found', len(eid_insertions), 'probe recordings from', brain_region)
```

found 43 probe recordings from LGD

Example06: Query Good-IBL units

```
In [2]: from one.api import ONE
        from brainbox.io.one import SpikeSortingLoader
        from ibllib.atlas import AllenAtlas

        one = ONE()
        ba = AllenAtlas()

        pid = 'da8dfec1-d265-44e8-84ce-6ae9c109b8bd'
        #LOAD SPIKES
        sl = SpikeSortingLoader(pid=pid, one=one, atlas=ba)
        spikes, clusters, channels = sl.load_spike_sorting()
        clusters = sl.merge_clusters(spikes, clusters, channels)

        #Filter GOOD - IBL - UNITS
        good_clusterIDs = clusters['cluster_id'][clusters['label'] == 1]
        good_clusterIDs
```


```
Out[2]: array([ 0,  1,  6, 35, 37, 38, 44, 49, 51, 53, 54, 63, 80,
                99, 109, 114, 127, 128, 129, 130, 132, 167, 169, 177, 188, 198,
                208, 217, 232, 243, 244, 249, 259, 261, 274, 288, 296, 297, 305,
                312, 313, 320, 325, 326, 327, 328, 339, 340, 341, 355, 356, 360,
                364, 367, 375, 378, 384, 393, 394, 395, 399, 401, 403, 406, 407,
                416, 425, 435, 436, 438, 446, 448, 452, 456, 458, 459, 461, 462,
                464, 469, 472, 475, 476, 479, 482, 488, 489, 494, 495, 500, 505,
                507, 509, 510, 518, 519, 527, 528, 529, 531, 536, 537, 543, 549,
                550, 554, 556, 560, 567, 570, 575, 577, 578, 581, 583, 587, 588,
                590, 594, 595, 597, 598, 599, 601, 607, 612, 613, 615, 619, 622,
                625, 626, 639, 641, 645, 648, 654, 665, 666, 667, 668, 669, 673,
                675, 677, 683, 685, 687, 715, 741, 743, 777, 782, 785, 792, 793,
                804, 812, 821, 826, 830, 856, 861, 869, 913], dtype=int64)
```

Criteria Good-IBL unit:


- 1) amplitude > 50 uV
- 2) noise cut-off < 20
- 3) refractory period violation

Practical examples


Example07_LoadWheel




Example08_LoadData_PlotPSTH



Example09_LoadData_PlotFR



Example10_ReactionTimes



IBL modular data architecture allows:

- organization of data within a lab
- integration of data from multiple labs
- Flexible access during project development
- pipelined analysis, display on an interactive website
- multiple access methods
- Modules can be used by individual labs, large or small collaborations

Useful links

ONE documentation

https://int-brain-lab.github.io/iblenv/notebooks_external/one_quickstart.html#

Viz website

<https://viz.internationalbrainlab.org/app>

Main website

<https://www.internationalbrainlab.com>

Task being performed by the mouse-

<https://doi.org/10.7554/eLife.63711>

Neural data that has been recorded

https://int-brain-lab.github.io/iblenv/notebooks_external/data_release_repro_ephys.html